



LW MEASUREMENTS LLC

Category	High Resolution Balances		
Model	HRB		
Model Suffix	103	203	303
Casing	Glass Draft Shield		
Capacity	100g	200g	300g
2nd Capacity	No		
Resolution	0.001g	0.001g	0.001g
No. of Display Divisions	100,000	200,000	300,000
No. of Internal Resolution	16m		
Accuracy	+/-5d		
Platter Size(inch)	3.85" Ø	3.85" Ø	3.85" Ø
Draft Shield	Glass		
Scale Dimensions(inch)	7.5 x 10 x 10.5		
LCD Height(inch)	1		
Beeper	Yes		
Auto Backlight	Yes,Selectable		
Auto Shut Off	Yes,Selectable		
Auto Hold	No		
Unit Switching	g,oz,lb,dwt,ozt,ct,tlt,tlh,tlj,gn,dr,mm,tola,gsm,T/A/R,T/M/R		
AD Processor	24 bit A/D		
AD Update Speed	7.5 Hz		
Settling Time	<2s		
Calibration Software	AnyCAL		
No. of Keys	5		
Memory	Zero		
Counting Function	No		
RS232C	Yes		
USB	No		
Battery	Yes, 6 x AA size battery		
AC Adapter	Yes, 12VDC		
Shipping Weight(lb)	11		
Pieces per master ctn	2		
Units per pallet	60		