



LW MEASUREMENTS LLC



Category	High Resolution Balance	
Model	HRB	
Model Suffix	10001LP	20001
Casing	Large Platter	
Capacity	10kg	20kg
2nd Capacity	No	
Resolution	0.1g	0.1g
No. of Display Divisions	100,000	200,000
No. of Internal Resolutions	16m	16m
Accuracy	+/-3d	
Platter Size(inch)	7.8 x 8.5	
Draft Shield	No	
Scale Dimensions(inch)	12.25 x 8.5 x 3.75	
LCD Height(inch)	1.1	
Beeper	Yes	
Auto Backlight	Yes, Selectable	
Auto Shut Off	Yes, Selectable	
Auto Hold	No	
Unit Switching	g,oz,lb,dwt,ozt,ct,tlt,tlh,tlj,gn,dr,mm,tola,gsm,T/A/R,T/M/R	
AD Processor	24 bit A/D	
AD Update Speed	7.5 Hz.	
Settling Time	<2s	
Calibration Software	AnyCal	
No of Keys	5	
Memory	Zero	
Counting Function	Yes	
RS232C	Yes	
USB	No	
Battery	Rechargeable	
AC Adapter	Yes, 10VDC	
Shipping Weight(lb)	10	
Pieces per master ctn	4	
Units per pallet	100	