



LW MEASUREMENTS LLC



Category	Magnetic Force Balances
Model	HRB
Model Suffix	4102
Casing	Plastic
Capacity	4100g
2nd Capacity	No
Resolution	10mg
No. of Display Divisions	410,000
No. of Internal Resolutions	10 million
Accuracy	e=100mg
Platter Size(inch)	185 x 185mm
Draft Shield	No
Scale Dimensions(inch)	12 x 8 x 3.5
LCD Height(inch)	0.56
Beeper	Yes
Auto Backlight	Yes, Selectable
Auto Shut Off	Yes, Selectable
Auto Hold	No
Unit Switching	
AD Processor	Yes
AD Update Speed	
Settling Time	<1.5s
Calibration Software	Single point CAL
No of Keys	7
Memory	Zero
Counting Function	Yes
RS232C	Yes
USB	Yes
Battery	No
AC Adapter	Yes, 7.5VDC
Shipping Weight(lb)	6.7
Pieces per master ctn	1
Units per pallet	60