

T SCALE

LW MEASUREMENTS LLC



Category	High Res Balance
Model	IHB
Model Suffix	PLUS 2
Casing	ABS plastic
Capacity	300g
2nd Capacity	No
Resolution	0.001g
No. of Display Divisions	300,000d
No. of Internal Resolutions	600,000
Accuracy	+/- 5d
Platter Size	3.15" round
Draft Shield	Glass
Scale Dimensions	13.4" x 8.7"x 12.6"
LCD Height	0.65"
Beeper	Yes
Auto Backlight	Yes
Auto Shut Off	Yes
Auto Hold	No
Unit Switching	g,ct,lb,oz,d,yn,ozt,dwt,mom,tlt,tlh,tlj,t,bt, n
AD Processor	24 bit A/D
AD Update Speed	15-60
Settling Time	<2S
Calibration Software	both lb and kg calibrations
No of Keys	7
Memory	Zero
Counting Function	Yes
RS232C	Yes
USB	No
Battery	Yes, 6V/1.2A
AC Adapter	Yes, 12V/500mA
Shipping Weight(lb)	10 lbs
Pieces per master ctn	4
Units per pallet	24 set