

	
LW MEASUREMENTS LLC	
<b>Category</b>	<b>Wheel Chair Scale</b>
<b>Model</b>	<b>LWC</b>
<b>Model Suffix</b>	800
<b>Casing</b>	Aluminum
<b>Capacity</b>	800lb
<b>2nd Capacity</b>	No
<b>Resolution</b>	0.2lb
<b>No. of Display Divisions</b>	4000
<b>No. of Internal Resolutions</b>	2m
<b>Accuracy</b>	+/-3d
<b>Platter Size(inch)</b>	37.8 x 39.4
<b>Draft Shield</b>	No
<b>Scale Dimensions(inch)</b>	43.9 x 39.4 x3
<b>LCD Height(inch)</b>	1
<b>Beeper</b>	Yes
<b>Auto Backlight</b>	Yes,Selectable
<b>Auto Shut Off</b>	Yes,Selectable
<b>Auto Hold</b>	No
<b>Unit Switching</b>	kg,lb
<b>AD Processor</b>	21 bit A/D
<b>AD Update Speed</b>	10 Hz
<b>Settling Time</b>	<2s
<b>Calibration Software</b>	AnyCAL
<b>No of Keys</b>	5
<b>Memory</b>	Zero
<b>Counting Function</b>	No
<b>RS232C</b>	Yes
<b>USB</b>	No
<b>Battery</b>	6 x AAsize battery
<b>AC Adapter</b>	Yes, 9VDC
<b>Shipping Weight(lb)</b>	78
<b>Pieces per master ctn</b>	1
<b>Units per pallet</b>	15