

LW MEASUREMENTS LLC



Category	Magnetic Force Balances	
Model	HRB	
Model Suffix	224	423
Casing	Glass Draft Shield	Glass Draft Shield
Capacity	220g	420g
2nd Capacity	No	No
Resolution	0.1mg	1mg
No. of Display Divisions	2.2m	420,000
No. of Internal Resolutions	10 million	10 million
Accuracy	e=1mg	e=10mg
Platter Size(inch)	85mm	110mm
Draft Shield	Glass	
Scale Dimensions(inch)	12 x 8 x 12	
LCD Height(inch)	0.56	
Beeper	Yes	
Auto Backlight	Yes, Selectable	
Auto Shut Off	Yes, Selectable	
Auto Hold	No	
Unit Switching	Yes	
AD Processor		
AD Update Speed		
Settling Time	<3s	<1.5s
Calibration Software	Single point CAL	
No of Keys	6	7
Memory	Zero	
Counting Function	Yes	
RS232C	Yes	
USB	Yes	
Battery	No	
AC Adapter	Yes, 7.5VDC	
Shipping Weight(lb)	7.8	
Pieces per master ctn	1	
Units per pallet		