



Category	K - Series Balances			
Model	KMR	KHR		
Model Suffix	6000	3000	500	123
Casing	Plastic			Plastic
Capacity	6000g	3000g	500g	120g
2nd Capacity	No	No	No	No
Resolution	1g	0.1g	0.01g	0.001g
No. of Display Divisions	6,000	30,000	50,000	120,000
No. of Internal Resolutions	2 m		2 m	4m
Accuracy	+/-3d		+/-3d	+/-3d
Platter Size(inch)	5.7 x 5.7		5.7 x 5.7	3.85" Ø
Draft Shield	Bowl	Bowl	Bowl	Mini Glass draft shield
Scale Dimensions(inch)	7 x 9 x 5		7 x 9 x 5	8.9 x 7.3 x 6.2
LCD Height(inch)	0.6		0.6	0.6
Beeper	Yes		Yes	Yes
Auto Backlight	Yes,Selectable		Yes,Selectable	Yes,Selectable
Auto Shut Off	Yes,Selectable		Yes,Selectable	Yes,Selectable
Auto Hold	No	No	No	No
Unit Switching	kg,g,lb,oz,lb:oz		g,lb,oz,lb:oz	g,mg,dwt,ozt,gn
AD Processor	21 bit A/D		21 bit A/D	22 bit A/D
AD Update Speed	7.5 Hz		10Hz	10Hz
Settling Time	<2s		<2s	<2s
Calibration Software	AnyCAL		AnyCAL	AnyCAL
No of Keys	4		4	4
Memory	Zero		Zero	Zero
Counting Function	No		No	No
RS232C	No		No	No
USB	No		No	No
Battery	6 x AA size battery		6 x AA size battery	6 x AA size battery
AC Adapter	Yes, 9V DC		Yes, 9V DC	Yes, 9V DC
Shipping Weight(lb)	3.2		3.2	4
Pieces per master ctn	6	6	6	4
Units per pallet	192	192	192	80