



LW MEASUREMENTS LLC



Category	High Resolution Balances	
Model	HRB	
Model Suffix	602	1002
Casing	Glass Draft Shield	
Capacity	600g	1000g
2nd Capacity	No	
Resolution	0.01g	0.01g
No. of Display Divisions	60,000	100,000
No. of Internal Resolutions	16m	
Accuracy	+/-3d	
Platter Size(inch)	4.7 x 5.3	
Draft Shield	Glass	
Scale Dimensions(inch)	7.5 x 10 x 10.5	
LCD Height(inch)	1	
Beeper	Yes	
Auto Backlight	Yes, Selectable	
Auto Shut Off	Yes, Selectable	
Auto Hold	No	
Unit Switching	g,oz,lb,dwt,ozt,ct,tlt,tlh,tlj,gn,dr,mm,tola,gsm,T/A/R,T/M/R	
AD Processor	24 bit A/D	
AD Update Speed	7.5 Hz	
Settling Time	<2s	
Calibration Software	AnyCAL	
No of Keys	5	
Memory	Zero	
Counting Function	No	
RS232C	Yes	
USB	No	
Battery	Yes, 6 x AA size battery	
AC Adapter	Yes, 12VDC	
Shipping Weight(lb)	11 with outer carton	
Pieces per master ctn	2	
Units per pallet	60	