

												
<b>Category</b>	<b>High Resolution Balances</b>											
<b>Model</b>	<b>HRB</b>											
<b>Model Suffix</b>	103	303	1002	1002TL	3002	3001	10001	10001LP	20001	30001	60001	
<b>Casing</b>	Glass Draft Shield			Plastic				Large Platter		Die Cast Aluminum		
<b>Capacity</b>	100g	300g	1000g	1000g	3000g	3000g	10kg	10kg	20kg	30kg	60kg	
<b>Resolution(kg)</b>	No			No				No		No		
<b>Resolution(g)</b>	0.001g		0.01g			0.1g			0.1g		0.2	
<b>Resolution(oz)</b>	0.00005	0.0001	0.0005		0.001		0.005			0.01		
<b>Resolution(lb)</b>	0.00001		0.00002			0.0002			0.0002			0-99.99lb/0.00
<b>Resolution(dwt)</b>	0.001		0.01			0.1			0.1		0.2	
<b>Resolution(ozt)</b>	0.00005		0.0005			0.005			0.005		0.01	
<b>Resolution(ct)</b>	0.005	0.01	0.05		0.1		0.5			1		
<b>Resolution(tl.t)</b>	0.00005		0.0005			0.005			0.005			
<b>Resolution(tl.h)</b>	0.00005		0.0005			0.005			0.005			
<b>Resolution(tl.J)</b>	0.00005		0.0005			0.005			0.005			
<b>Resolution(gn)</b>	0.02		0.2			2			2		5	
<b>Resolution(dr)</b>	0.001		0.01		0.1		0.1			0.2		
<b>Resolution(mm)</b>	0.0005		0.005			0.05			0.05			
<b>Resolution(tola)</b>	0.0001		0.001		0.01			0.01		0.02		
<b>Resolution(gsm)</b>	0.1		1			10			No			
<b>Resolution(tar)</b>	0.1		1			1			No			
<b>Resolution(tma)</b>	0.1		1			1			No			
<b>Display Divisions</b>	100,000	300,000	100,000	300,000	30,000	100,000	200,000	300,000				
<b>Internal Resolutions</b>	16m											
<b>Accuracy</b>	+/-5d		+/-3d		+/-5d		+/-3d		+/-3d		+/-3~5d	
<b>Platter Size(inch)</b>	3.85" Ø		3.85" Ø		4.7 x 5.3		6.5 x 7		7.8 x 8.5		12.4 x 12.7	
<b>Draft Shield</b>	Glass			No				No				
<b>Scale Dimensions(inch)</b>	7.5 x 10 x 10.5			7.5 x 10 x 3.5				12.25 x 8.5 x 3.75		12.4 x 12.7 x 3.2		
<b>LCD Height(inch)</b>	1			1				1.1		1		
<b>Beeper</b>	Yes											
<b>Auto Backlight</b>	Yes,Selectable											
<b>Auto Shut Off</b>	Yes,Selectable											
<b>Auto Hold</b>	No											
<b>Unit Switching</b>	g,oz,lb,dwt,ozt,ct,tlh,tlj,gn,dr,mm,tola,gsm,T/A/R,T/M/R											
<b>AD Processor</b>	24 bit A/D											
<b>AD Update Speed</b>	10 hz											
<b>Settling Time</b>	<2s											
<b>Calibration Software</b>	AnyCAL											
<b>No of Keys</b>	5											
<b>Memory</b>	Zero											
<b>Counting Function</b>	No			Yes				Yes				
<b>RS232C</b>	Yes											
<b>USB</b>	No											
<b>Battery</b>	Yes, 6 x AA size battery						Rechargeable					
<b>AC Adapter</b>	Yes, 12VDC											
<b>Shipping Weight(lb)</b>	11 with outer carton			6				10		18.5		
<b>Pieces per master ctn</b>	2			6				4		2		
<b>Units per pallet</b>	60			144				100		84		