



LW MEASUREMENTS



Category	Vet Scales		
Model	LVS / LVS+	LVS-XL / LVS-XL+	LVS-S
Model Suffix	700		
Casing	Steel		
Capacity	700lb		
Resolution(kg)	0.1		
Resolution(lb)	0.2		
Resolution(oz)	No		
Display Divisions	3500		
Internal Resolutions	2 m		
Accuracy	+/-3d		
Platter Size(inch)	37.4 x 19.7	44 x 22	37.4X19.7
Scale Dimensions(inch)	37.4 x 19.7 x 3	44 x 22 x 3	37.4X19.7X2.8
LCD Height(inch)	1		
Beeper	Yes		
Auto Backlight	Yes,Selectable		
Auto Shut Off	Yes,Selectable		
Auto Hold	No		
Unit Switching	kg,lb		
AD Processor	21 bit A/D		
AD Update Speed	10 Hz		
Settling Time	<2s		
Calibration Software	AnyCAL		
No of Keys	5		
Memory	Zero		
Counting Function	No		
RS232C	Yes		
USB	No		
Battery	6 x AAsize battery		
AC Adapter	Yes, 9VDC		
Shipping Weight(lb)	47	58	99
Pieces per master ctn	1		
Units per pallet	20	20	15