



Category	LOW COST VET SCALE		
Model	LC-VS 60	LC-VS 330	LC-VS 400
Model Suffix	60	330	400
Casing	Tempered Glass		Plastic
Capacity	60lb	330lb	400lb
2nd Capacity	No		
Resolution(kg)	0.01	0.05	
Resolution(lb)	0.02	0.1	
No. of Display Divisions	3000	3300	4000
No. of Internal Resolutions	2 million		
Accuracy	+/- 3~5 d	+/- 3d	
Platter Size(inch)	25 x 16	33.1 x 16	37.4 x 19.7
Scale Dimensions(inch)	25 x 16 x 1.9	33.1 x 16 x 1.9	37.4 x 19.7 x 3
LCD Height(inch)	1.1		1
Beeper	Yes		
Auto Backlight	Yes, Selectable		
Auto Shut Off	Yes, Selectable		
Auto Hold	No		
Unit Switching	kg, lb., lb. Oz	kg, lb., lb. Oz	kg, lb.
AD Processor	21 bit A/D		
AD Update Speed	10 Hz		
Settling Time	<2s		
Calibration Software	Not AnyCal		AnyCal
Applicable calibration weights	lbs or kg		
No of Keys	4		5
Memory	Zero		
Counting Function	No		
RS232C	Yes		
USB	No		
Battery	4 x AA Battery	4 x AA battery	6 x AA Battery
AC Adapter	6V/300mah	6V/300mah	9V/500mah
Shipping Weight(lb)	17lbs	23lbs	36lbs
Operating environment	0°C - 40°C / 32°F - 104°F (no condensation)		
Standard accessories	Manual, AC adapter, rubber protective mat		Manual, AC adapter, rubber protective mat, indicator
Accessories sold separately	Wired indicator (sold separately)		None