



Category	Jewelry Balance		Mid Resolution Balances				
Model	JLY53	JLY123	MRB				
Model Suffix	53	123	500	1200	2500	6000	10000
Casing	Top Loader Balances						
Capacity	50g	120g	500g	1200g	2500	6000g	10,000g
Resolution(kg)	NO		0.0001		0.001		
Resolution(g)	0.001g		0.1		1		
Resolution(oz)	0.00005		0.005		0.05		
Resolution(lb)	No		0.0002		0.002		
Resolution(dwt)	0.001		No				
Resolution(ozt)	0.00005		No				
Resolution(ct)	0.005		No				
Resolution(tl.t)	No						
Resolution(tl.h)	No						
Resolution(tl.J)	No						
Resolution(gn)	No						
Resolution(dr)	No						
Resolution(mm)	0.0005		No				
Resolution(tola)	No						
Resolution(gsm)	No						
Resolution(tar)	No						
Resolution(tma)	No						
Display Divisions	50,000	120,000	5000	12000	25000	6000	10000
Internal Resolutions	16m		2 m				
Accuracy	+/-3d	+/-5d	+/-3d				
Platter Size(inch)	2.67" Ø		6.25 x 7				
Draft Shield	No						
Scale Dimensions(inch)	6.1 x 4.5 x 3.5		0.6				
LCD Height(inch)	0.4		0.6				
Beeper	YES						
Auto Backlight	Yes,Selectable						
Auto Shut Off	Yes,Selectable						
Auto Hold	No						
Unit Switching			g,oz,lb				
AD Processor	24 bit A/D		21 bit A/D				
AD Update Speed	7.5 Hz		10 Hz				
Settling Time			<2s				
Calibration Software	Single point CAL		AnyCAL				
No of Keys	3		4				
Memory			Zero				
Counting Function	YES		No				
RS232C	No						
USB	No						
Battery	yes,6XAAA size battery		Yes, 6 x AA size battery				
AC Adapter	Yes, 9VDC						
Shipping Weight(lb)	2		3				
Pieces per master ctn	8		12				
Units per pallet	384		288				