



Category	Wheel Chair Scale		
Model	LWC 1000	LWC-HR	LWC-P
Model Suffix	1000		800
Casing	Aluminum	Steel	Aluminum
Capacity	1000 lb		800 lb
Resolution(kg)	0.1		
Resolution(lb)	0.2		
Resolution(oz)	No		
Display Divisions	5000		4000
Internal Resolutions	2m		
Accuracy	+/-3d		
Base Size(inch)	39.4 X 37	39.4 x 37	37.25 x 6.25
Platter size(inch)	36.5 x 31	36.5 x 31	37.25 x 6.25
Scale Dimensions(inch)	43.9 x 39.4 x 3	45 x 39.4 x 3.1	54.3 x 7.5 x 3
LCD Height(inch)	1		
Beeper	Yes		
Auto Backlight	Yes,Selectable		
Auto Shut Off	Yes,Selectable		
Auto Hold	No		
Unit Switching	kg,lb		
AD Processor	21 bit A/D		
AD Update Speed	10 Hz		
Settling Time	<2s		
Calibration Software	AnyCAL		
No of Keys	5		
Memory	Zero		
Counting Function	No		
RS232C	Yes		
USB	No		Yes
Battery	6 x AA size battery		74 V rechargeable Li
AC Adapter	Yes, 9 VDC		Yes, 12 VDC
Shipping Weight(lb)	100	140	43 (Each Rail 16 lb)
Pieces per master ctn	1		
Units per pallet	10		28