



LC-VS Series

Low-Cost Vet Scales



LC-VS 60



Compact!

LC-VS 330



Versatile!

LC-VS 400



Accurate!

LW Measurements LLC takes pride in producing quality products that are affordable, durable, and most importantly accurate. Visit our website for more information and videos of our products.

www.lwmeasurements.com



The Low-Cost Vet Scale or LC-VS

Are you looking for a cost-effective solution to accurately measure the weight of animals in your veterinary practice? Look no further! Our low-cost vet scales are designed to provide reliable and accurate weight measurements for animals of all sizes. Whether you're a small clinic or a larger practice, our scales are the perfect addition to your veterinary equipment lineup.

Key Features

Precision Accuracy: Our vet scales are engineered to provide precise weight measurements, ensuring you have accurate data for treatment and record-keeping.

Durable Build: Built to withstand the demands of a busy veterinary practice, our scales are constructed with high-quality materials to ensure durability and longevity.

Various Weight Capacities: Choose from a range of weight capacities to accommodate animals of different sizes, from small pets to larger livestock.

Easy-to-Read Display: The large, easy-to-read display ensures that weight readings can be quickly and accurately assessed.

Tare Functionality: The tare function allows you to subtract the weight of blankets or containers, giving you the true weight of the animal.

Tempered glass for easy cleaning (LC-VS 60/330)

Less effort cleaning and sturdy platform for all animals.



Buttons for pushing with fingers and toes...

The button design on the LC-VS 60 and 330 were meant for fingers at counter height and the toe of your sneaker if you preferred to set it up at ground level. This was designed to bring convenience to you and your decision of where to place this scale.



Convenience, portability, and capability are important.

Each scale is provided with AC power Adapters, rubber protective mats, adjustable feet, and wheels (LC-VS 400). The LC-VS 400 comes with a Remo indicator that is factory set for the specific capacity of this scale model. The LC-VS 60/330 has the capability of plugging in a remote indicator (sold separately) which is manufactured by LW Measurements. All three scales models offer portability. They run off batteries as well as a RS232C (LC-VS 400 only) cable port for connecting to computers and printers.

Compact and Low Profile

The LC-VS 60 and 330 models are low profile and small enough to set on a counter or big enough to lay on the floor. Features AC power source or batteries and a port to plug in a remote indicator.



AC adapter port



Remote Indicator connection port

LCVS 60 and 330 remote indicator (sold separately)



Caster wheels




Versatile and Portable

Adjustable footing

Companion Animals: Our scales are perfect for weighing small to medium-sized pets, ensuring accurate dosing of medications, and monitoring of their health.

Livestock: From farm animals to equine patients, our vet scales provide an efficient and accurate way to monitor the weight of larger animals.

Exotic Species: Our scales are versatile enough to handle a wide range of animals, including exotics. Provide the best care for all your patients.



Stainless Steel Carrying Handle

The Remo Indicator...

Allows for easy readability and functionality. Hold, Zero, Tare, and Print at the press of a button. Even the ability to connect a printer or to a computer.

Features

LC-VS 60

- 60lb capacity
- Auto Shut off
- Tempered glass base
- Backlight function
- AnyCal calibration
- Zero function
- Units: lbs / kg (option)
- Rubber protective mat for base
- Easy to read LCD display
- AC adapter or 4 x AA batteries

LC-VS 330

- 330lb capacity
- Auto Shut off
- AnyCal calibration
- Tempered glass base
- Rubber protective mat
- Zero function
- Units: lbs / kg (option)
- Backlight function
- Easy to read LCD display
- AC adapter or 4 x AA batteries

LC-VS 400

- 400lb capacity
- Durable plastic base
- Steel frame
- REMO indicator
- Auto shut off
- AnyCal calibration
- Auto backlight function
- Units: lbs / kg (option)
- Stain-less handle & wheels
- AC adapter or 6 x AA batteries

Specifications

LW Measurements LLC	LC-VS 60	LC-VS 330	LC-VS 400
Model Suffix	60	330	400
Casing	Tempered Glass		Plastic
Capacity	60lb	330lb	400lb
2nd Capacity	No		
Resolution(kg)	0.01	0.05	
Resolution(lb)	0.02	0.1	
No. of Display Divisions	3000	3300	4000
No. of Internal Resolutions	2 million		
Accuracy	+/- 3~5 d	+/- 3d	
Platter Size(inch)	25 x 16	33.1 x 16	37.4 x 19.7
Scale Dimensions(inch)	25 x 16 x 1.9	33.1 x 16 x 1.9	37.4 x 19.7 x 3
LCD Height(inch)	1.1		1
Beeper	Yes		
Auto Backlight	Yes, Selectable		
Auto Shut Off	Yes, Selectable		
Auto Hold	No		
Unit Switching		kg, lb., lb. Oz	kg, lb.
AD Processor	21 bit A/D		
AD Update Speed	10 Hz		
Settling Time	<2s		
Calibration Software	Not AnyCal		AnyCal
Applicable calibration weight	lbs or kg		
No of Keys	4		5
Memory	Zero		
Counting Function	No		
RS232C	Yes		
USB	No		
Battery	4 x AA Battery	4 x AA battery	6 x AA Battery
AC Adapter		6V/300mah	9V/500mah
Shipping Weight(lb)	17lbs	23lbs	36lbs
Operating environment	0°C - 40°C / 32°F - 104°F (no condensation)		
Standard accessories	Manual, AC adapter, rubber protective mat		Manual, AC adapter, rubber protective mat, indicator
Accessories sold separately	Wired indicator (sold separately)		None

Warranty

LW Measurements products are under warranty against factory defects for a period of two (2) years from the date of shipment for Customers in the continental United States.

LW Measurements will pay for freight both ways within the first 30 days of purchase if technical support is required. (Determined by our certified Technical Support Team) In the event technical support is needed the customer is responsible for shipping the product back to us after the 30-day period has expired. Once the product is received, the technical support team will inspect and, as necessary, repair or replace the scale. LW Measurements will ship the product back to the customer at our expense.

Any new scales returned for warranty must be properly packaged in the original box. If they are not properly packed in the original box, the customer is responsible for the shipping cost. If we determine there is a factory defect, we will provide a replacement and accept responsibility for shipping back to the customer. If we determine that it is not a factory defect, the customer will be responsible for shipping fees incurred.

For Customers outside the lower 48 States, including Mexico, Canada, Puerto Rico, Hawaii, Alaska, and all other countries are responsible for all shipping fees incurred.

Our warranty does not cover misuse or neglect including but not limited to battery or water damage, overloading, and chewed or cut wires, keypad damage (punctures), scratched platforms, dented platforms unless caused by shipping (Must report/return within 30 days). If the product is found to have been misused or damaged by the customer, LW Measurements is not responsible for the cost of return.

For warranty claims please call your customer service representative.



LW Measurements LLC

Visit www.lwmeasurements.com or call us at 1-707-542-2185

