



<b>Category</b>	<b>Analytical Balance</b>	
<b>Model</b>	<b>DAB / ES-J</b>	
<b>Model Suffix</b>	224 / 220	224A / 220A
<b>Casing</b>		
<b>Capacity</b>	220g	220g
<b>Resolution(Kg)</b>		
<b>Resolution(g)</b>	0.0001g	
<b>Resolution(oz)</b>	0.00000353	
<b>Resolution(ct)</b>	0.0005	
<b>Resolution(lb)</b>	0.00000022	
<b>Resolution(ozt)</b>	-	-
<b>Resolution(dwt)</b>	-	-
<b>Resolution(tL.t)</b>	-	-
<b>Resolution(tL.h)</b>	-	-
<b>Resolution(tL.J)</b>	-	-
<b>Resolution(gn)</b>	-	-
<b>Resolution(dr)</b>	-	-
<b>Resolution(mm)</b>	-	-
<b>Resolution(tola)</b>	-	-
<b>Resolution(gsm)</b>	-	-
<b>Resolution(tar)</b>	-	-
<b>Resolution(tma)</b>	-	-
<b>Display Divisions</b>	2200000	
<b>Internal Resolutions</b>	4	
<b>Accuracy</b>	+/- 3d	
<b>Platter Size(inch)</b>	3.5 " Ø	
<b>Draft Shield</b>		
<b>Scale Dimensions(inch)</b>		
<b>LCD Height(inch)</b>	1.15	
<b>Beeper</b>		
<b>Auto Backlight</b>		
<b>Auto Shut Off</b>		
<b>Auto Hold</b>		
<b>Unit Switching</b>	g, lb, oz, ct	
<b>AD Processor</b>		
<b>AD Update Speed</b>		
<b>Settling Time</b>		
<b>Calibration Software</b>	Auto Internal Calibration	
<b>No of Keys</b>	6	
<b>Memory</b>		
<b>Counting Function</b>		
<b>RS232C</b>		
<b>USB</b>		
<b>Battery</b>		
<b>AC Adapter</b>	Yes, 12V DC	
<b>Shipping Weight(lb)</b>		
<b>Pieces per master ctn</b>		
<b>Units per pallet</b>	24	