

					
<b>Category</b>	<b>High Precision Balance</b>				
<b>Model</b>	<b>DPB / AP</b>		<b>DPB / CL</b>	<b>DPB / BP</b>	
<b>Model Suffix</b>	323 / 3203	623 / 6203	502 CL / 1002	32002	62002
<b>Casing</b>	Glass Draft Shield			Plastic	
<b>Capacity</b>	320g	620g	1000g	3200g	6200g
<b>Resolution(Kg)</b>	No				
<b>Resolution(g)</b>	0.001g		0.01g	0.01g	
<b>Resolution(oz)</b>	0.0000353		0.000353	0.000353	
<b>Resolution(ct)</b>	0.005		0.05	0.05	
<b>Resolution(lb)</b>	-	-	-	-	-
<b>Resolution(ozt)</b>	-	-	-	-	-
<b>Resolution(dwt)</b>	-	-	-	-	-
<b>Resolution(tl.t)</b>	-	-	-	-	-
<b>Resolution(tl.h)</b>	-	-	-	-	-
<b>Resolution(tl.J)</b>	-	-	-	-	-
<b>Resolution(gn)</b>	-	-	-	-	-
<b>Resolution(dr)</b>	-	-	-	-	-
<b>Resolution(mm)</b>	-	-	-	-	-
<b>Resolution(tola)</b>	-	-	-	-	-
<b>Resolution(gsm)</b>	-	-	-	-	-
<b>Resolution(tar)</b>	-	-	-	-	-
<b>Resolution(tma)</b>	-	-	-	-	-
<b>Display Divisions</b>	320000	620000	1000000	320000	620000
<b>Internal Resolutions</b>	3				
<b>Accuracy</b>	+/- 5d				
<b>Platter Size(Inch)</b>	4.5 " Ø		3.5 " Ø	6.3 " Ø	
<b>Draft Shield</b>	Glass			No	
<b>Scale Dimensions(Inch)</b>	8.27x 12.2 x 12.4		7.32 x 10.83 x 10.43	8.27x 12.2 x 3.54	
<b>LCD Height(Inch)</b>	1.34		0.86		
<b>Beeper</b>	No				
<b>Auto Backlight</b>	Yes , Selectable				
<b>Auto Shut Off</b>	No				
<b>Auto Hold</b>	No				
<b>Unit Switching</b>	g, oz,ct				
<b>AD Processor</b>	No				
<b>AD Update Speed</b>	No				
<b>Settling Time</b>	No				
<b>Calibration Software</b>	No				
<b>No of Keys</b>	8		6	8	
<b>Memory</b>	Zero				
<b>Counting Function</b>	Yes				
<b>RS232C</b>	Yes				
<b>USB</b>	Yes				
<b>Battery</b>	No				
<b>AC Adapter</b>	Yes, 7.5V DC		Yes, 12V DC	Yes, 7.5V DC	
<b>Shipping Weight(lb)</b>	15			13	
<b>Pieces per master ctn</b>	1			4	
<b>Units per pallet</b>	30		72	64	